

**Baxter**

# Acute RRT Modalities: Comparisons and Considerations





## Table of Contents

**1**

**ACUTE KIDNEY INJURY (AKI) OVERVIEW**

**2**

**RRT MODALITIES FOR AKI**

**3**

**RRT MODALITY CONSIDERATIONS**

**4**


**SUMMARY**

**5**

**ACRONYMS/ABBREVIATIONS/REFERENCES**



Acute kidney injury is **common** among hospitalized patients globally<sup>1</sup>



AKI AFFECTS AN ESTIMATED **20%** OF HOSPITALIZED PATIENTS WORLDWIDE<sup>1,\*</sup>

AKI is a **serious** condition


AKI IS ASSOCIATED WITH AN INCREASED RISK OF **MORBIDITY and MORTALITY**<sup>2-6</sup>

AKI IS ASSOCIATED WITH AN INCREASED RISK OF **CKD, including ESRD**<sup>7-9</sup>


\*Multicentre meta-analysis of 154 studies (n=3,585,911), primarily in hospital settings, that adopted a KDIGO-equivalent AKI definition between 2004 and 2012. Pooled rates.<sup>1</sup>




**Fluid overload** is one condition that may **adversely** impact AKI patient prognosis<sup>10,11</sup>



Consequences of fluid overload may lead to organ dysfunction<sup>13</sup>



Fluid overload at RRT initiation for AKI has been associated with an increased risk of mortality<sup>11,\*</sup>



**FLUID OVERLOAD** IN PATIENTS WITH **AKI** IS A **SERIOUS CONDITION**<sup>14-16</sup>


\*Prospective, observational cohort study of 296 adults treated with RRT in 17 Finnish ICUs from Sep 2011–Feb 2012.<sup>11</sup>



## AKI is associated with **substantial financial burden**<sup>17-19,\*</sup>

### AKI status impacts daily costs<sup>18</sup>

PATIENTS WITH AKI HAVE SIGNIFICANTLY HIGHER DAILY COSTS COMPARED WITH PATIENTS WITHOUT AKI<sup>18,†</sup>



### AKI is expensive even relative to other acute medical conditions<sup>19</sup>

Acute medical condition	Adjusted mean cost difference, in 2012 USD (95% CI) <sup>a</sup>
<b>AKI-D<sup>b</sup></b>	<b>11,016 (10,468, 11,564)</b>
Sepsis	4822 (4696, 5068)
VTE	3782 (3611, 3953)
Acute pancreatitis	1802 (1676, 1929)
<b>AKI<sup>c</sup></b>	<b>1795 (1692, 1899)</b>
Pneumonia	1705 (1584, 1825)
Stroke	1427 (1281, 1573)
MI	14 (-91, 119)
GI bleed	-860 (-961, -759)

THE INCREMENTAL COST OF AKI-D OR AKI IS HIGHER THAN FOR MANY OTHER CONDITIONS FOUND IN HOSPITALIZED PATIENTS<sup>19,‡</sup>

<sup>a</sup>Compared with reference group without the condition of interest.

<sup>b</sup>Compared with patients without AKI. <sup>c</sup>Includes patients with dialysis-requiring AKI (AKI-D).

WHILE EXPENDITURES MAY VARY BY COUNTRY,  
**AKI is a COSTLY CONDITION**<sup>17-19</sup>


\*Costs for hospitalisation due to AKI may vary from country to country.

†Multicentre, retrospective cohort study of 659,945 adult hospital admissions across central China in 2013.<sup>18</sup>

‡2012 multicentre, retrospective study of 29,763,649 adult US hospitalisations without ESRD.<sup>19</sup>



Various **renal replacement modalities** are available for the **management of AKI**<sup>20-24</sup>



ACUTE RRT IS DELIVERED AS **EITHER**  
A **CONTINUOUS OR**  
**INTERMITTENT** THERAPY<sup>20</sup>



Modalities **differ** in their typical characteristics<sup>25</sup>

Typical RRT modality characteristics and settings for a 70-kg AKI patient<sup>25-27</sup>

Parameter	CONTINUOUS THERAPIES			INTERMITTENT THERAPIES	
	CVVH	CVVHD	CVVHDF	SLED*	IHD
Blood flow (Q <sub>B</sub> , mL/min)	150-250	150-250	150-250	100-300	200-300
Predominant solute transport principle					
Ultrafiltrate (mL/h)	1500-2000	variable	1000-1500	variable	variable
Dialysate flow (Q <sub>D</sub> , mL/h)	0	1500-2000	1000-1500	6000-18,000	18,000-30,000
Replacement fluid for zero balance (mL/h)	1500-2000	0	1000-1500	0	0
Urea clearance (mL/min)	25-33	25-33	25-33	80-90	200-500

\*SLED is a type of PIRRT.<sup>21</sup>



Convection




Diffusion

Q<sub>B</sub>, Q<sub>D</sub>, AND UREA CLEARANCE TEND TO BE **LOWER** IN **CONTINUOUS** THERAPIES THAN IN **INTERMITTENT** THERAPIES<sup>25-27</sup>



**Individual patient needs** can be addressed by considering the characteristics of the various **RRT modalities**<sup>28</sup>

Relative **features, risks, and burdens** of different RRT modalities<sup>28</sup>



EACH RRT MODALITY HAS POTENTIAL **BENEFITS** AND **LIMITATIONS** FOR THE MANAGEMENT OF PATIENTS WITH AKI<sup>28</sup>





**Selection of RRT modality** requires careful consideration of many patient- and ICU-specific factors<sup>25,28</sup>

## Overview of modality considerations



**CLINICAL CONSIDERATIONS:  
FLUID OVERLOAD AND HEMODYNAMIC INSTABILITY**



**CLINICAL CONSIDERATIONS: LONG-TERM OUTCOMES**



**MACHINE AND PRESCRIPTION CONSIDERATIONS**



**SOLUTION CONSIDERATIONS**



**LONG-TERM COST CONSIDERATIONS**



**EQUIPMENT FOOTPRINT AND MOBILITY CONSIDERATIONS**



## Clinical considerations: fluid overload and hemodynamic instability




Fluid overload in AKI patients can be treated by fluid removal during RRT, but rapid fluid removal that does not allow time for plasma refill may lead to **hemodynamic instability**<sup>25,29</sup>



**Avoiding rapid fluid removal** to prevent hypovolaemia may **improve** AKI patient outcomes<sup>25,29</sup>

## Modality comparisons<sup>30</sup>



**CRRT IS A PREFERRED RRT**  
 BY MANY CLINICIANS FOR AKI PATIENTS WHO ARE  
**HEMODYNAMICALLY UNSTABLE**<sup>25,29</sup>




## Clinical considerations: long-term outcomes



AKI is associated with an **increased risk** of long-term dialysis dependence;<sup>8</sup> acute RRT modality type may impact this risk<sup>31-34</sup>

Patients on chronic dialysis at day 90 by initial RRT modality<sup>31,\*</sup>



CHRONIC DIALYSIS HAZARD RATIO (95% CI) FOR CRRT vs IHD WAS 0.75 (0.65-0.87), P<0.0001<sup>31</sup>

## Modality comparisons<sup>31-34</sup>

CONTINUOUS	INTERMITTENT	
CRRT	PIRRT	IHD
Patients are <b>less likely to require chronic dialysis</b> following initial AKI episode compared with patients treated with IHD	Insufficient evidence	It has been reported that patients are <b>more likely to require chronic dialysis</b> following initial AKI episode compared with patients treated with CRRT

USE OF **CRRT** FOR AKI MANAGEMENT HAS BEEN ASSOCIATED WITH A **LOWER RISK of CHRONIC DIALYSIS** COMPARED WITH IHD<sup>31-34</sup>

\*Retrospective multicentre cohort study of critically ill adults with AKI between 1996 and 2009. 2004 patients originally treated with CRRT and 2004 patients originally treated with IHD were propensity matched and rates of dialysis dependence were compared.<sup>31</sup>




## Machine and prescription considerations



RRT machines deliver **different** dose intensities over **different** durations of therapy<sup>21,22,25</sup>

### Kinetic modeling of urea clearance by different RRT modalities<sup>35</sup>



A SAWTOOTH PATTERN WAS OBSERVED WHEN USING INTERMITTENT THERAPIES TO REMOVE UREA, WHILE CONTINUOUS THERAPY MAINTAINED A CONSISTENT BUN LEVEL OVER TIME<sup>35</sup>

### Modality comparisons<sup>21,22,35</sup>

CONTINUOUS	INTERMITTENT	
CVVH	SLED	IHD
<p>Intended to run 24 h/day</p> <ul style="list-style-type: none"> <li>Slow but continuous urea clearance helps avoid spikes in BUN levels</li> </ul>	<p>Typically run in 6–12 h sessions delivered 3–7 times/week</p> <ul style="list-style-type: none"> <li>Intermittent nature does not allow for continuous urea clearance, which could result in variable BUN levels</li> </ul>	<p>Typically run in 3–6 h sessions delivered 3–7 times/week</p> <ul style="list-style-type: none"> <li>Intermittent nature does not allow for continuous urea clearance, which could result in variable BUN levels</li> </ul>

UNLIKE IHD OR PIRRT, **CRRT** IS RUN ON MACHINES THAT DELIVER **CONTINUOUS** SOLUTE REMOVAL<sup>22,35</sup>



## Solution considerations



Typically, CRRT solutions are **commercially** prepared, while IHD and PIRRT use **local water sources** to prepare dialysate<sup>29,36,37</sup>



Preparing solutions on-line from local water sources **necessitates** water **treatment** and routine water **quality monitoring** to assure clean water standards are met<sup>36-38</sup>

## Modality comparisons<sup>29,36-40</sup>

CONTINUOUS	INTERMITTENT	
CRRT	PIRRT	IHD
<p>Because no on-line solutions are typically used, <b>no water treatment systems</b> are required</p> <ul style="list-style-type: none"> <li>Monitoring water quality is not applicable</li> </ul>	<p>If a centralized water treatment system is unavailable in the ICU, <b>individual water quality monitoring is necessary</b></p> <ul style="list-style-type: none"> <li>If a centralized water treatment system is not used, <b>staff need to monitor dialysate quality</b> for individual patients</li> <li><b>Disinfection requirements</b> may limit treatment duration to &lt;12 hours<sup>41</sup></li> </ul>	<p>If a centralized water treatment system is unavailable in the ICU, <b>individual water quality monitoring is necessary</b></p> <ul style="list-style-type: none"> <li>If a centralized water treatment system is not used, <b>staff need to monitor dialysate quality</b> for individual patients</li> </ul>

WATER TREATMENT AND QUALITY TESTING MAY CONTRIBUTE TO **INCREASED MONITORING** WHEN USING SOLUTIONS PREPARED ON-LINE FOR **IHD and PIRRT**<sup>39,42</sup>




## Long-term cost considerations



Because **initial RRT modality** may impact the risk of **chronic dialysis**,<sup>31</sup> **long-term costs** of AKI may also be influenced by **initial treatment modality**<sup>43</sup>


Cumulative costs of dialysis dependence by initial AKI treatment modality<sup>43,\*</sup>



**MEAN 5-YEAR TOTAL COST/PATIENT OF AKI-D<sup>+</sup> WAS \$37,780 FOR CRRT AS THE INITIAL MODALITY COMPARED WITH \$39,448 FOR IRRT<sup>43</sup>**

\*Including cost of dialysis dependence. Cost in 2013 USD.

## Modality comparisons<sup>43</sup>



THE **LONG-TERM COST** OF AKI MAY BE **LOWER** FOR PATIENTS INITIALLY TREATED WITH **CRRT** COMPARED TO THOSE TREATED WITH IHD<sup>43</sup>

\*Health outcomes and healthcare costs were simulated and averaged for a cohort of 1000 patients initiated on CRRT and a cohort of 1000 patients initiated on IRRT. All costs were inflated to 2013 USD.<sup>43</sup>



## Equipment footprint and mobility considerations



Water treatment systems required for IHD and PIRRT add to physical **space** requirements and water lines may limit RRT mobility in ICUs without central water treatment systems<sup>37,40</sup>



In ICUs without central water treatment systems, portable **water treatment devices** may be necessary,<sup>40</sup> which can occupy as much as 0.13–0.16 m<sup>2</sup> of floor space<sup>44,45</sup>

## Modality comparisons<sup>36,37,40,41,46–49</sup>

CONTINUOUS	INTERMITTENT	
CRRT	PIRRT	IHD
<p>Because the CRRT machine is the <b>only component</b> that contributes to the therapy's physical footprint, treatment mobility may be increased</p> <ul style="list-style-type: none"> <li>•</li> </ul> <p>No space considerations for water treatment systems are necessary</p>	<p>Both the IHD machine and water treatment systems contribute to the therapy's physical footprint, which may <b>impact treatment mobility</b> in ICUs without central water treatment systems</p> <ul style="list-style-type: none"> <li>•</li> </ul> <p>In situations where a central water treatment system is not utilised, the <b>greater physical footprint</b> of the machine + water treatment system may impact ICU spacing</p>	

**WATER TREATMENT** EQUIPMENT MAY ADD TO THE FOOTPRINT OF **IHD** AND **PIRRT** SYSTEMS, POTENTIALLY DECREASING TREATMENT MOBILITY AND IMPACTING SPACING CONSIDERATIONS<sup>40,47–49</sup>



AKI is a **common** and **costly** condition among ICU patients,<sup>1,17-19</sup> and is associated with increased risks of **morbidity and mortality**<sup>2-9</sup>



Acute RRT is delivered as **either a continuous or an intermittent** therapy, each of which have unique characteristics, settings, and limitations<sup>20,25-28</sup>

**Selection of RRT modality** requires careful consideration of many patient- and ICU-specific factors<sup>25,28</sup>



FLUID OVERLOAD AND  
HEMODYNAMIC  
INSTABILITY



LONG-TERM  
CLINICAL  
OUTCOMES



MACHINE AND  
PRESCRIPTION



SOLUTIONS



LONG-TERM COSTS



EQUIPMENT  
FOOTPRINT AND  
MOBILITY

**CRRT** IS A PREFERRED  
RENAL REPLACEMENT THERAPY  
BY MANY CLINICIANS FOR PATIENTS WITH AKI  
WHO ARE HEMODYNAMICALLY UNSTABLE<sup>25,28</sup>







AKI, acute kidney injury; AKI-D, dialysis-requiring AKI; BUN, blood urea nitrogen; CI, confidence interval; CKD, chronic kidney disease; CRRT, continuous renal replacement therapy; CVVH, continuous veno-venous hemofiltration; CVVHD, continuous veno-venous hemodialysis; CVVHDF, continuous veno-venous hemodiafiltration; dL, decilitre; EDHF, extended daily hemofiltration; ESRD, end-stage renal disease; Feb, February; GI, gastrointestinal; h, hour; HD, hemodialysis; ICU, intensive care unit; IHD, intermittent hemodialysis; IRRRT, intermittent renal replacement therapy; K<sup>+</sup>, potassium ion; KDIGO, Kidney Disease Improving Global Outcomes; kg, kilogram; m<sup>2</sup>, square meters; mg, milligram; MI, myocardial infarction; min, minute; mL, millilitre; PD, peritoneal dialysis; PIRRT, prolonged intermittent renal replacement therapy; NH<sub>3</sub>, ammonia; Q<sub>B</sub>, blood flow rate; Q<sub>D</sub>, dialysis flow rate; RRT, renal replacement therapy; Sep, September; SLED, sustained or slow low-efficiency dialysis; SLEDD, sustained or slow low-efficiency daily dialysis; US, United States; USD, United States dollar; vs, versus; VTE, venous thromboembolism

1. Susantitaphong P, et al. *Clin J Am Soc Nephrol*. 2013;8:1482-1493.
2. Coca SG, et al. *Am J Kidney Dis*. 2009;53:961-973.
3. Ricci Z, et al. *Kidney Int*. 2008;73:538-546.
4. Chawla LS, et al. *Clin J Am Soc Nephrol*. 2014;9:448-456.
5. Brown JR, et al. *Ann Thorac Surg*. 2016;102:1482-1489.
6. Wu VC, et al. *J Am Soc Nephrol*. 2014;25:595-605.
7. Ishani A, et al. *J Am Soc Nephrol*. 2009;20:223-228.
8. Coca SG, et al. *Kidney Int*. 2012;81:442-448.
9. Wald R, et al. *JAMA*. 2009;302:1179-1185.
10. Bellomo R, et al. *Crit Care Med*. 2012;40:1753-1760.
11. Vaara ST, et al. *Crit Care*. 2012;16:R197.
12. Kaddourah A, et al. *N Engl J Med*. 2017;376:11-20.
13. Prowle JR, et al. *Nat Rev Nephrol*. 2010;6:107-115.
14. Zhang L, et al. *J Crit Care*. 2015;30:860.e7-13.
15. Bouchard J, et al. *Kidney Int*. 2009;76:422-427.
16. Heung M, et al. *Nephrol Dial Transplant*. 2012;27:956-961.
17. National Clinical Guideline Centre (UK). Acute Kidney Injury: Prevention, Detection and Management Up to the Point of Renal Replacement Therapy. 2013; NICE Clinical Guidelines, No. 169. Introduction.
18. Xu X, et al. *Clin J Am Soc Nephrol*. 2015;10:1510-1518.
19. Silver SA, et al. *J Hosp Med*. 2017;12:70-76.
20. Fleming GM. *Organogenesis*. 2011;7:2-12.
21. O'Reilly P, Tolwani A. *Crit Care Clin*. 2005;367-378.
22. Pannu N, Gibney RTN. *Ther Clin Risk Manag*. 2005;1:141-150.
23. Sun Z, et al. *Crit Care*. 2014;18:R70.
24. Kitchlu A, et al. *BMC Nephrol*. 2015;16:127.
25. Kidney Disease: Improving Global Outcomes (KDIGO) Acute Kidney Injury Work Group. KDIGO Clinical Practice Guideline for Acute Kidney Injury. *Kidney Int Suppl*. 2012;2:1-138.
26. Baxter Corporation. REVACLEAR Dialyzer Technology. 2017. Available from: [http://www.baxter.ca/en\\_CA/assets/downloads/2017/Revaclear%20Spec%20Sheet%20Brochure%20English.pdf](http://www.baxter.ca/en_CA/assets/downloads/2017/Revaclear%20Spec%20Sheet%20Brochure%20English.pdf) (accessed December 2018).
27. Fresenius Medical Care. Optiflux Dialyzers. 2016. Available from: [http://www.fmcna-dialyzers.com/images/pdf/101046-Optiflux-Dialyzer\\_SpecSheet.pdf](http://www.fmcna-dialyzers.com/images/pdf/101046-Optiflux-Dialyzer_SpecSheet.pdf) (accessed December 2018).
28. Ostermann M, et al. *Blood Purif*. 2016;42:224-237.
29. Murugan R, et al. *Blood Purif*. 2016;42:266-278.
30. Ostermann M. In: Bellomo R, et al (eds). 40 Years of Continuous Renal Replacement Therapy. *Contrib Nephrol*. 2018;194:51-59.
31. Wald R, et al. *Crit Care Med*. 2014;42:868-877.
32. Bell M, et al. *Intensive Care Med*. 2007;33:773-780.
33. Cartin-Ceba R, et al. *Intensive Care Med*. 2009;35:2087-2095.
34. Lin YF, et al. *Am J Surg*. 2009;198:325-332.
35. Liao Z, et al. *Artif Organs*. 2003;27:802-807.
36. Coulliette AD, Arduino MJ. *Semin Dial*. 2013;26:427-438.
37. Glorieux G, et al. *Nephrol Dial Transplant*. 2012;27:4010-4021.
38. Azar AT, Ahmad S. Hemodialysis Water Treatment System. In: Azar A (eds). Modelling and Control of Dialysis Systems. Studies in Computational Intelligence. 2013;404:347-378.
39. BC Renal Agency. Clinical Practice Standards and Procedures for Dialysis Water Quality: 2b: Endotoxin Testing of Dialysis Water. 2011. Available from: [http://www.bcrenalagency.ca/resource-gallery/Documents/2bEndotoxin-Testing-of-Dialysis-Water-Final\\_2012.pdf](http://www.bcrenalagency.ca/resource-gallery/Documents/2bEndotoxin-Testing-of-Dialysis-Water-Final_2012.pdf) (accessed December 2018).
40. Bellomo R, et al. *Crit Care Resusc*. 2002;4:281-290.
41. Fresenius Medical Care. 2008K<sup>2</sup> Hemodialysis Machine Operator's Manual. 2016. Available from: [https://fmcna.com/wp-content/uploads/documents/490136\\_Rev\\_J.pdf](https://fmcna.com/wp-content/uploads/documents/490136_Rev_J.pdf) (accessed December 2018).
42. Food and Drug Administration. Quality Assurance Guidelines for Hemodialysis Devices. 1991. Available from: <https://www.fda.gov/downloads/medicaldevices/deviceregulationandguidance/guidancedocuments/ucm073435.pdf> (accessed December 2018).
43. Ethgen O, et al. *Nephrol Dial Transplant*. 2015;90:54-61.
44. Mar Cor Purification. 700 Series Portable. 2014. Available from: [http://www.mcpur.com/main/library/12\\_brochures/1238642\\_\(700\).pdf](http://www.mcpur.com/main/library/12_brochures/1238642_(700).pdf) (accessed December 2018).
45. Mar Cor Purification. Millenium HX. 2015. Available from: [http://www.mcpur.com/main/library/12\\_brochures/3027573\\_\(MHX\).pdf](http://www.mcpur.com/main/library/12_brochures/3027573_(MHX).pdf) (accessed December 2018).
46. Baxter Healthcare Corporation. The PRISMAFLEX System. 2016. Available from: [https://www.baxter.com/sites/g/files/ebysai746/files/2017-11/Prismaflex-07.11-Brochure-New\\_Accts.pdf](https://www.baxter.com/sites/g/files/ebysai746/files/2017-11/Prismaflex-07.11-Brochure-New_Accts.pdf) (accessed December 2018).
47. Ledebro I, Blankestijn PJ. *NDT Plus*. 2010;3:8-16.
48. Kasperek T, Rodriguez OE. *Clin J Am Soc Nephrol*. 2015;10:1061-1071.
49. Poeppl K, et al. *Vet Clin Small Anim*. 2011;177-191.